Generalidades

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General information

Read these instructions carefully before using the product. Failure to comply with the safety and storage instructions may cause harm to human health and to the environment.

Intended use

The MYU P/S/F activator is a medical device supplied in non-sterile packaging. It acts gently on the structure of the oral cavity, improving the functions that take place within it: swallowing, breathing, chewing and phonation.

Warnings

The activator has been prescribed by your dentist; therefore, we advise you to follow all instructions carefully to obtain the best result in the shortest time possible.

- Before use, make sure that there are no specific allergies to the material. In the extremely rare event that allergic reactions occur, stop using it and contact your dentist.
- Biting down too hard on the activator may cause it to break, either partially or completely. If this happens, stop using it and contact your dentist.
- When disposing of the activator and its packaging, strictly comply with the applicable waste disposal regulations.
- Keep the product and the manual out of the reach of children and pets.

In the very unlikely event that a serious incident occurs in relation to the use of this device, it must be reported to the manufacturer, who will in turn notify the competent authority of the Member State to which the user is subject.

This instruction manual must be kept throughout the use of the product.

Methods of use

The dentist chooses the activator model that is most suitable for the patient, and only the dentist may modify it to make it more functional for the therapy.

Remove the activator from the case in which it is to be stored and moisten it before inserting it between the dental arches. If it is the first use, rinse it with water or clean it more thoroughly with specific natural cleaning products. Then check that the insertion position is correct.

If, during use, you experience problems such as pain, ulcers, inflammation or swelling of the gums, stop using it and contact your dentist.

Dosage

The dentist will determine the wearing times. Normally, the activator should be used at night while sleeping — during the first few days it is normal to lose it during sleep or to wake up with tooth pain — and during the day for at least 2 hours.

During daytime use:

- **The patient must not play or talk with the activator in the mouth.**
- **The lips must be kept closed and the tip of the tongue in contact with the palate.**
- **For mouth-breathing patients, it will be necessary to practise breathing through the nose.**
- For patients with atypical swallowing, it will be necessary to practise swallowing correctly while keeping the activator in the mouth.
- **Exercises with the activator may be performed, individualised for each patient and indicated by the dentist.**

Age of use

The ideal age for achieving the best results with the activator is between 6 and 15 years. From the age of 16 onwards, excellent results can still be obtained, but usually over longer periods of time. The activator may be used from the time the deciduous dentition develops, and there is no age limit for its

use. It is also recommended for use as retention instead of a retainer, or after the use of fixed orthodontic appliances.

Contraindications

Carefully consider the use of the activator in cases of mobile teeth, advanced periodontal disease and pain in the temporomandibular joint, until the dentist has assessed suitability for use.

Maintenance

The quality of the material used and the smooth shapes make cleaning easy with water or, more thoroughly, with specific natural cleaning products. Once clean, make sure that both the activator and its case remain dry.

Replacement criteria

It is up to the dentist to decide when to replace the activator. In general, this occurs when a pre-established therapeutic objective has been achieved. Replacement is recommended when there are structural changes due to chewing that may compromise its function. In any case, the activator must be replaced within a maximum period of 12 months of use.

Storage

Keep the activator in a dry, clean place, away from sources of heat or direct sunlight; do not place it in contact with citrus fruits or essential oils, as these may stain it.

Always make sure that the hygienic conditions in which the activator is stored are appropriate for its use.

Avoid introducing dirt and germs into your mouth through the appliance: always handle it with clean hands and store it properly.

Materials used

The activator is manufactured from thermoplastic material composed of the following components:

- Hydrogenated styrene block copolymer: SEBS (CAS 66070-58-4), SEPS/SEEPS (CAS 132778-07-5)
- Thermoplastic olefin polymers: PP (CAS 9003-07-0, 9010-79-1), PE (CAS 9002-88-4), EVA (CAS 24937-78-8)
- Styrenic thermoplastic polymers: PS (CAS 9003-53-6), SBC (CAS 9003-55-8)
- Paraffinic plasticising oils: CAS 8042-47-5, 8012-95-1
- Mineral filler: calcium carbonate (CAS 471-34-1)
- Various additives

If desired, you may request the material safety data sheet from your dentist.

Expiry date

The device, if stored in its packaging and, in any case, in an environment suitable for preservation, may be used within a period of 5 years.

Explanation of symbols

SYMBOL	NAME	DESCRIPTION
	Manufacturer	Indicates the manufacturer of the medical device as stipulated in European Directives 90/385/EEC, 93/42/EEC and 98/79/EC. This symbol is accompanied by the name and address of the manufacturer. Este símbolo se acompaña del nombre y la dirección del fabricante.
	European Conformity	European Conformity (CE) mark for Class I medical devices.
	Catalogue number	Indicates the manufacturer's catalogue number used to identify the medical device.
	Date of manufacture	Indicates the date of manufacture of the medical device. Date format: MM/YYYY Formato de fecha: MM/AAAA
	Batch code	Indicates the manufacturer's batch code used to identify the batch or series.

	Non-sterile	Indicates that the medical device has not undergone a sterilisation process.
	Consult instructions for use	Indicates that the user must consult the instructions for use.
	Medical device	Indicates that the product is a medical device.

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